

MMC HF

High-frequency surgery Units

High Class Devices • Modern Technology (MOSFET)

MMC 100 S

- With coagulation mode
- 110W power
- For small surgeries
- Easy and safe use
- Monopolar and bipolar Mode



MMC 200S

- For more comfort and more important surgeries
- 180W power
- Easy and safe use
- More efficiency in coagulation
- Monopolar and bipolar Mode

MMC 400S

- The tool of the surgeon
- 360W, digital control
- Double switch handle for cutting / coagulation
- Independent regulation of the level of hemostasis
- Automatic contact plate control
- Simultaneous monopolar and bipolar connection
- Advanced technology



MMC surgery units

Fields of application from high-frequency surgery in veterinary medicine:

- Hemostasis
- Epulidectomy
- Gingivectomy
- Tumor ablation
- Mammectomy
- Surgeries in the oral cavity
- Orthopedic surgeries...

Benefits of the high-frequency surgery:

- Marginal bleedings because of coagulation effect
- Good healing process
- Profitable, less need of suture material
- High success rate from the surgeries
- Large fields of application

Technical features:

MMC 100S

Dimensions: 90 x 225 x 160 mm
3.5'' x 9'' x 6''

Weight: 2.2 kg / 5 lb

Line voltage: 110/120 V /
220/240 V

Line frequency: AC 50/60 Hz

Monopolar output:

- Cut: 110W
- Coagulation: 80W

Bipolar output:

- Cut: 100W
- Coagulation: 75W

Frequency: 300 KHz

MMC 200S

Dimensions: 90 x 225 x 160 mm
3.5'' x 9'' x 6''

Weight: 2.2 kg / 5 lb

Line voltage: 110/120 V /
220/240 V

Line frequency: AC 50/60 Hz

Monopolar output:

- Cut: 180W
- Coagulation: 158W

Bipolar output:

- Cut: 162W
- Coagulation: 141W

Frequency: 300 KHz

MMC 400S

Dimensions: 130x330x280 mm
5'' x 13'' x 11''

Weight: 5.2 kg / 11.5 lb

Line voltage: 110/120 V /
220/240 V

Line frequency: AC 50/60 Hz

Monopolar output:

- Cut: 360W
- Coagulation: 200W

Bipolar output:

- Cut: 90W
- Coagulation: 50W

Frequency: 300 KHz

Accessories:

- 4 electrodes (needle, blade, ball and loop)
- Cables
- Double foot switch
- Neutral electrode (optional rubber plate)
- Bipolar forceps (optional)



Medical equipment for veterinarians

15, rue du Noroit
22100 Taden - France

Tél. +33 (0)2 96 85 86 76

Fax +33 (0)2 96 85 86 77

welcome@manomedical.com

www.manomedical.com