




# SK-500I V2 Syringe Pump

For an ACCURATE  
and SAFE infusion



**Mano  
Médical**

Medical equipment for veterinarians

# SK-500I V2

## Syringe pump



### Application:

- During anaesthesia
- In analgesia
- In oncology



### Flow and performance :

- Compatible with 5, 10, 20, 30, 50 and 60 ml syringes
- VTBI: 0,1 - 999,9 ml
- Flow rate range:
  - ▶ 5 ml syringe: 0.1 - 150 ml / h
  - ▶ 10 ml syringe: 0.1 - 300 ml / h
  - ▶ 20 ml syringe: 0.1 - 600 ml / h
  - ▶ 30 ml syringe: 0.1 - 900 ml / h
  - ▶ 50/60 ml syringe: 0.1 - 1500 ml / h
- Accuracy: +/- 2%
- Bolus rate: 0,10 - 1500 ml/h
- KVO rate (Keep Vein Open): 0.5 ml/h
- Occlusion:
  - ▶ Low: 300 ±125 mmHg
  - ▶ Medium: 525 ±125 mmHg
  - ▶ High: 900 ±180 mmHg



### Display:

- Flow rate
- Accumulated injection volume
- Preset volume
- Syringe size
- Battery capacity
- AC power indicator
- Battery charge indicator



### Alarm:

- VTBI done, VTBI near done
- KVO finish
- Occlusion
- Syringe empty, syringe near empty, syringe not installed
- System error, system abnormal
- Battery low, battery empty
- AC power disconnection



Dimensions:	360 x 145 x 123 mm
Weight:	1.95 Kg
Screen:	2.5'' LCD screen
Syringe sizes:	5, 10, 20, 30, 50/60 ml
Flow rate:	1 - 1500 ml / h
Flow rate limit :	1 - 9999 ml
Accuracy:	+/- 2%
Battery:	Lithium
Battery operation time:	> 6 hours
Charging time:	< 5 hours @100%
Power:	100-240 V, 50/60 Hz



Medical equipment for veterinarians

15, rue Noroit  
22100 Taden - France  
Tél. +33 (0)2 96 85 86 76  
Fax +33 (0)2 96 85 86 77

[welcome@manomedical.com](mailto:welcome@manomedical.com)

[www.manomedical.com](http://www.manomedical.com)