

DENTISTRY

Veterinary Specific Dentistry Devices

▶ Strong Micromotor including polishing function

Type E micromotor - light, compact and handy for interventions on small animals or devitalized teeth



- Regular speed from 0 to 35 000 turns/min.
- 1:1 contra angle
- Setting drive
- Compact and light type E motor
- Including handpiece troy and pedal



▶ Vet 200 Plus Ultrasonic scaler

Piezo-electric scaler - portable and powerful, perfect for veterinary use!



- Titanium base alloy handpiece
- 5 piezo-ceramic tips
- Automatic frequency tracking digitally controlled
- Easy installation
- Compact and light design
- CE certified

Including tips:



 **Mano Médical**
Medical equipment for veterinarians

15, rue du Noroit
22100 Taden - France
Tél. +33 (0)2 96 85 86 76

welcome@manomedical.com

www.manomedical.com