

Poskom X-ray Generator

Battery powered • LED collimator • Dual laser

- ▶ High frequency generator
- ▶ Internal battery or external power supply operation
- ▶ Up to 650 exposures by just one charging
- ▶ Dual laser for easy focal distance setting
- ▶ Strong body against external shock and contamination
- ▶ Fully charged only within 4 hours



Output	Power	90 kV and 20 mA
	mAs	0,4 - 20 mAs
	kV	50 -90 kV
	Max. output	1,6 kW @ 80 kV
X-ray tube	Focal spot	1.2 mm
Collimator	Focal distance	35 x 35 cm@ cm SID
	Lamp	LED
	Laser pointer	Built-in dual laser
Operation panel	Memory storage	15 APR
Battery	Charging time	4 hours
Weight		7.2 Kg
Dimensions		203 x 174 x 307 mm

 **Mano Médical**
Medical equipment for veterinarians

15, rue du Noroit
22100 Taden - France
Tél. +33 (0)2 96 85 86 76

welcome@manomedical.com

www.manomedical.com