

The new state-of-the-art MicroFlow™ capnograph specifically developed for veterinarians!

# MicroFlow™ Capnograph

## CAPNO-5



# MICROFLOW™ SIDE-STREAM CAPNOGRAPHY



Developed for point of care capnography monitoring of intubated or non-intubated patients in surgery, dentistry, procedural sedation, pre-op preparation, anesthesia recovery, and intensive / critical care, the Capno-5 provides:

- Real-time CO<sub>2</sub> waveforms (capnogram)
- End-tidal CO<sub>2</sub> (ETCO<sub>2</sub>)
- Fractional inspired CO<sub>2</sub> (FiCO<sub>2</sub>)
- Respiratory rate (RR)
- Audible and visual apnea alarms
- Delivered with table mount



With the very low sample flow rate of 50 ml per minute, the Capno-5 is suitable for use in a wide range of veterinary patients, including cats and dogs of all sizes!

And unlike all mainstream capnographs, which are not designed for use in dentistry and non-intubated patients, the Capno-5 can be used in all procedures for intubated or non-intubated patients.

## Technical features:

The Capno5 includes:

- Moisture trap
- Airway Adapter
- AC power adapter
- Table mount

Data ranges

- Sample flow rate: 50 ml / min.
- Respiration range: 2 -150 bpm
- Dimensions: 155 x 68 x 35 mm
- Weight: 340 gr
- Display: 3.5'' (8.9 cm)
- Battery autonomy: up to 10 hours



- ✓ EASY-TO-USE  
CAPNOGRAPH
- ✓ LONG-LIFE BATTERY  
(UP TO 10 HOURS OF USE)
- ✓ LARGE DISPLAY  
3.5'' (8.9 CM)



Mano  
Médical

Medical equipment for veterinarians

15, rue du Noroît  
22100 Taden - France  
Tel. +33 (0)2 96 85 86 76

welcome@manomedical.com

www.manomedical.com