


HF1060 X-ray generator

Mobile and powerful X-ray generator

Perfect both for your X-ray room and mobile use

- ▶ High frequency generator
- ▶ Exposure hand switch
- ▶ Ergonomic and soft operation panel layout
- ▶ Dual laser spot
- ▶ Convenient handling for advanced image quality
- ▶ Supplied in a robust carrying case



Output	Adjustment	kV and mAs
	mAs	0,32 - 100 mAs
	Max. output power	3,2 kW @ 80 kV 40 mA
	Accuracy	kV: +/-5%, mAs: +/-10%
X-ray tube	Focal spot	2.0 mm
Collimator	Focal spot	66 cm
Generator	Output	60 mA @ 40-60 kV 30 mA @ 81-100 kV
HF inverter		High Frequency
Operation panel	Display	1) kV, mAs 2) Ready, X-ray, wait, error LED
	Operation switch	1) kV and mAs 2) Dual Laser & Lamp ON/OFF 3) Display reverse 4) APR memory storage 5) 8 memories (Up to 16 memories)
Weight		12.9 Kg




Mano Médical
 Medical equipment for veterinarians

15, rue du Noroit
 22100 Taden - France
 Tel. +33 (0)2 96 85 86 76

welcome@manomedical.com