

• Oxygen Concentrator

VET OXYFLOW OCI

HIGH
PRESSURE
TO SUIT MOST OF
THE ANAESTHETIC
MACHINES!




Mano
Médical
by génia

Veterinary Oxygen Concentrator



The output pressure superior to most of other devices on the market allows to provide oxygen efficiently to all types of anaesthesia machine.



VetOxyFlow OC1 has been developed to meet the needs of veterinary clinics that do not have any piping system, and/or to be directly connected to a ventilator or anaesthetic machine. This compact generator requires very little space and delivers 5 L/min oxygen with exceptional performance.



- ✓ **Fully automatic**
- ✓ **Easy operation**
- ✓ **Compact and silent**
- ✓ **Easy maintenance**

Technical data:

▶ Purity oxygen	≥ 93% ±3% (V/V)
▶ Outlet pressure	1.2 - 1.4 Bar
▶ Flow rate	0-5 L/min
▶ Sound level	< 45 dB
▶ Power consumption	300W
▶ Voltage	AC220V±22V, 50 Hz±1Hz AC110V±15V, 60 Hz±1Hz
▶ Net weight	22 kg
▶ Dimensions	49 x27 x56 cm



15, rue du Noroît
22100 Taden - France
Tel. +33 (0)2 96 85 86 76

welcome@manomedical.com

www.manomedical.com