

# petMAP™

Blood Pressure Measurement Device

A family of veterinary blood pressure devices tailored for your specific clinical needs



Advanced oscillometric BP measurement device designed for use in the broadest range of patient sizes and blood pressures



Mano Médical  
by génia

# Choose the right model to suit your needs!

	PetMAP™ graphic II	PetMAP™ g3	PetMAP™+ II	PetMAP™ + II SpO <sub>2</sub>
Measured parameters	Blood pressure (systolic, diastolic and MAP)  Heart rate	Blood pressure (systolic, diastolic and MAP)  Heart rate  Temperature (esophageal probe sold separately)	Blood pressure (systolic, diastolic and MAP)  Heart rate  Temperature  ECG (optional module)  Capnography (optional module)	Blood pressure (systolic, diastolic and MAP)  Heart rate  Temperature  Pulse oxymetry  ECG (optional module)  Capnography (optional module)
Display	4.3'' (10.9 cm) touchscreen	5'' (12.7 cm) touchscreen		
Cuffs inflation	Automatic	Automatic and BP cycle time setting for automatically repeated BP measurements		
Power	4 AA Alkaline	4 AA NiMH and AC/DC adapter		
Included accessories	7 cuffs from 2 to 5.5 cm		7 cuffs from 2 to 5.5 cm  Esophageal probe for temperature	7 cuffs from 2 to 5.5 cm  Esophageal probe for temperature  SpO <sub>2</sub> lingual and toe sensors

## The options:

- Large size cuffs: 6.5 - 8 - 10 and 13 cm
- Protection cover
- Table and pole mounts

More information about our petMAP



Table mount



Pole mount



Protective cover



Mano  
Médical  
by génia

15, rue du Noroît  
22100 Taden - France  
Tel. +33 (0)2 96 85 86 76

welcome@manomedical.com

www.manomedical.com