

QuickCuff™
*Patent Pending

**INNOVATIVE BLOOD
PRESSURE CUFFS**

**MEASURING
BLOOD
PRESSURE HAS
NEVER BEEN
SO EASY!**




**mano
médical**
EQUIPMENT FOR VETS
by génia

Click and Go Easy Placement!

A revolutionary, easy to use and stress free way to measure blood pressure


- ✓ No velcro
- ✓ One-handed application
- ✓ Resistant
- ✓ Easy to clean
- ✓ Fits all different limbs (tapered & untapered)
- ✓ 1 color for each size (12 sizes from 2 to 13 cm)
- ✓ Optional storage rack



1 Select a QuickCuff
Color-coded for easy identification
12 QuickCuffs to choose from size 2.0 cm to 13 cm



2 Squeeze to tighten around limb
Sizing tabs secure cuff in place
If it fits snug tight, it's the right size



3 Squeeze again to release



Quick Cuff™ size	Circumference	
	Min	Max
2 cm	3.14 cm	4.40 cm
2.5 cm	4.08 cm	5.65 cm
3 cm	5.34 cm	6.28 cm
3.5 cm	5.97 cm	7.23 cm
4 cm	6.60 cm	8.80 cm
4.5 cm	7.54 cm	9.42 cm
5 cm	7.85 cm	10.37 cm
5.5 cm	8.80 cm	11 cm
6.5 cm	11 cm	13.19 cm
8 cm	14.45 cm	17.28 cm
10 cm	21.05 cm	23.88 cm
13 cm	24.82 cm	29.53 cm



petMAP™
Blood Pressure Measurement Device

Still available!

Have a look at the petMAP™ range on
www.manomedical.com/en/diagnostic



22 rue Charles de Gaulle
22100 Lanvallay - France
Tel. +33 (0)2 96 85 86 76

export@manomedical.com

www.manomedical.com