

Radiolucent Table

MOBILE • ROBUST • CONVENIENT

OPTIONAL
STAND FOR
X-RAY UNIT -
LARGE MOVABLE
DRAWER



Veterinary X-ray table



- ✓ A lubrication-free sliding system allows the drawer and stand to be easily moved along the length of the table.
- ✓ The table is mounted on casters with brakes, which enable its easy mobility.
- ✓ Under the table, a movable drawer covered with non-slip mat is attached. Its large size is suitable for 43x43 cm flat panel at least.
- ✓ An anti-slip mat is also available for the table top.



Convenient:

- ▶ Space-saving and greater comfort with the sliding drawer and X-ray unit stand: no need anymore to handle the animal on the table or move the table top
- ▶ Very manoeuvrable thanks to its 4 lockable casters
- ▶ Height adjustment of the arm to set the distance between the X-ray portable unit and the flat panel
- ▶ Extremely quiet operation

Robust:

- ▶ Any type of portable X-ray unit up to 35 kg
- ▶ Steel base, covered with powder coated paint resistant to chipping and scratching
- ▶ Maintenance-free and insensitive to the effects of dirt, water, chemicals, heat or impacts

Easy to maintain:

- ▶ Easy to keep clean and sterile
- ▶ Waterproof, full resistance to moisture and water vapor



Compatible with our generators:

- ▶ Orange HF1060
- ▶ Poskom 20BT Lite

Table dimensions: 120 (L) x 63 (l) x 80 (H) cm

Internal dimensions of the drawer: 47 x 57 cm

Adjustable X-ray unit stand : From 60 to 125 cm from the drawer surface

Anti-slip mat: Included for both drawer and table top

X-ray unit max. weight: 35 kg

Animal max. weight: 150 kg

HPL radiolucent table top thickness: 16 mm

Table weight: 56 kg without X-ray unit stand
70 kg with X-ray unit stand



22 rue Charles de Gaulle
22100 Lanvallay - France
Tel. +33 (0)2 96 85 86 76

contact@manomedical.com

www.manomedical.com