

Infusion Range



Fluidotherapy by Mano Médical

ManoPump iV1 Plus*

Infusion pump



Large touchscreen • Optional fluid warmer

- High-performance battery: up to 10 hours run time
- Large touchscreen: 4.3'' (10.75 cm)
- 5 occlusion levels (50, 150, 300, 600 and 900 mmHg)
- Up to 2 000 entries in the library
- 5 000 review data
- IV stand support and kennel mount included
- Optional fluid warmer module

ManoPump iV2

Infusion pump



Simple and intuitive control panel

- Up to 5 hours operating time
- 2.8'' (7 cm) color TFT screen
- Front control panel
- 4 occlusion levels (150, 300, 600 and 900 mmHg)
- 30 entries in the library
- 5 000 review data
- IV stand support and kennel mount included

ManoPump SP1

Syringe pump



High-performance battery • Standard syringes compatible

- Automatic recognition of syringe sizes (5, 10, 20, 30 & 50/60 ml)
- High-performance battery: up to 12 hours run time
- Large touchscreen: 4.3'' (10.75 cm)
- Modes: ml/h (time, rate mode), body-weight and micro
- Flow rate: 1 500 ml/h
- 5 000 review data
- IV stand support included

Vet Temp3

Fluid warmer



3 lengths of heating grooves for flexibility and accuracy

- 3 heating groove lengths: 25, 31 and 52 cm
- Exclusive IV set or standard PVC sets of 20 and 60 drops/ml
- 3.4 and 4.5 mm IV set outer diameter
- Overheat & low temperatures, system error, door open, heating reminder alarms
- IV stand lanyard and bracket included

- [More information on
www.manomedical.com/en/infusion](http://www.manomedical.com/en/infusion)

*ManoPump iV 1 model also available
without fluid warmer module option



22 rue Charles de Gaulle
22100 Lanvallay - France
Tel. +33 (0)2 96 85 86 76

contact@manomedical.com

www.manomedical.com

Choose the model that suits your needs!

	ManoPump iV1	ManoPump iV1 Plus	ManoPump iV2
	Ref. POM200	Ref. POM210	Ref. POM300
	Up to 10 hours run time on battery		Up to 5 hours run time on battery
	Display on a 4.3'' (10.75 cm) screen		Display on a 2.8'' (7 cm) screen
	Touchscreen		Control buttons under display
	5 occlusion levels (50-900 mmHg)		4 occlusion levels (150-900 mmHg)
	Up to 2 000 entries in the library		30 entries in the library
	5 000 review data		2 000 review data
	Flow rate from 0.1 to 2 000 ml/h		
	IV stand support and kennel mount		
	-	Optional fluid warmer module (ref. POM202)	-
Infusion Pumps			

	ManoPump SP1		Vet Temp3
	Ref. SYR200		Ref. WAR163
	Up to 12 hours run time on battery		3 heating groove lengths : 25, 31 and 52 cm
	Display on a 4.3'' (10.75 cm) screen		Compatible with sets of 20 and 60 drops/ml
	Touscreen		3.4 and 4.5 mm IV set outer diameters
	Automatic recognition of syringe sizes (5, 10, 20, 30 and 50/60 ml)		Overheat & low temperatures, system error, door open and heating reminder alarms
	Modes: ml/h (time, rate mode), body-weight and micro		IV stand lanyard and bracket attaching methods
	5 000 review data		
	Flow rate from 0.1 to 1 500 ml/h		
	IV stand support and kennel mount		
Syringe pump		Fluid warmer	

22 rue Charles de Gaulle
22100 LANVALLAY
France

Tel. +33 (0)2 96 85 86 76
export@manomedical.com



www.manomedical.com

